

# Great White<sup>®</sup> Platinum Stirrups

*Gynecology, laparoscopy  
& urology procedures  
are safer than ever*

The Great White<sup>®</sup> Platinum Stirrups provide safe & quick positioning assistance for improved surgical site access while offering exceptionally precise placement.

## Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -30° low to 75° high; Adduction -9° and abduction +25°
- The raised lateral fin boots effectively minimize the possibility of nerve damage
- Included Platinum Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

**#800-0342-PL** Great White<sup>®</sup> Platinum Stirrups

### **Required Accessory: Purchase Two SpringLoc Clamp**

**#800-0338**

**#800-0338-EU**

**#800-0338-JPN**

**#800-0338-SWISS**



**600 lb. (272 kg)  
Patient Weight  
Capacity**



Patent  
Pending

### **Optional Accessories**

**#800-0074-S** Stirrup Dolly

**#800-0284** Stirrup Wall Rack

### **Replacement Pads**

**#508-1415** Platinum Stirrup Boot Pads, Set

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351