

Great White® Premium Stirrups

High-quality, feature-rich & priced right!

Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning

Features

- Exclusive 3-axis control from handle mimics the anatomy of the human hip
- Lithotomy Range -30° low to 75° high; Adduction -9° and abduction +25°
- Included premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-PR Premium Stirrups
#800-0342-PRF Premium Stirrups w/Fins

Required Accessory: Purchase Two SpringLoc Clamp

#800-0338 **#800-0338-EU**
#800-0338-JPN **#800-0338-SWISS**

Optional Accessories

#800-0074-S Stirrup Dolly
#800-0284 Stirrup Wall Rack
#800-0160 Clear View Leg Drapes, Sterile, 10/cs*
 *Sterile products are non-returnable.

Replacement Pads

#508-1354 Premium Stirrup Boot Pads, Set
#508-1415 Platinum Stirrup Boot Pads, Set



Patent
 Pending



#800-0342-PRF