

Great White® Robotic Stirrups

Robotic Gynecologic & Urologic procedures have never been safer

Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -55° low to 60° high; Adduction -9° and abduction +25°
- Included premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-R Great White® Robotic Stirrups
#800-0342-RF Great White® Robotic Stirrups w/Fins

Required Accessory: Purchase Two SpringLoc Clamp

#800-0338 **#800-0338-EU**
#800-0338-JPN **#800-0338-SWISS**

Optional Accessories

#800-0284 Stirrup Wall Rack
#800-0074-R Robotic Stirrup Dolly
#800-0160 Clear View Leg Drapes, Sterile, 10/cs*

*Sterile products are non-returnable.


Replacement Pads

#508-1354 Premium Stirrup Boot Pads, Set
#508-1415 Platinum Stirrup Boot Pads, Set



Patent Pending

600 lb. (272 kg)
Patient Weight Capacity



#800-0342-RF

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351