

# Great White<sup>®</sup> Robotic Stirrups

*Robotic Gynecologic & Urologic procedures have never been safer*

## Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -55° low to 60° high; Adduction -9° and abduction +25°
- Included premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- Lithotomy indicators help to accurately reposition legs

**#800-0342-R** Great White<sup>®</sup> Robotic Stirrups  
**#800-0342-RF** Great White<sup>®</sup> Robotic Stirrups w/Fins

## Required Accessory: Purchase Two SpringLoc Clamp

**#800-0338**                      **#800-0338-EU**  
**#800-0338-JPN**              **#800-0338-SWISS**

## Optional Accessories


**#800-0284** Stirrup Wall Rack  
**#800-0074-R** Robotic Stirrup Dolly

## Replacement Pads

**#508-1354** Premium Stirrup Boot Pads, Set  
**#508-1415** Platinum Stirrup Boot Pads, Set



**600 lb. (272 kg)**  
**Patient Weight**  
**Capacity**



#800-0342-RF

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351