

Great White® Robotic Stirrups

Robotic Gynecologic & Urologic procedures have never been safer

Features

- Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- Lithotomy Range -60° low to 40° high; Adduction -9° and abduction +25°
- Included Premium Boot Pads wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- Lithotomy indicators help to accurately reposition legs

#800-0342-RF Great White® Robotic Stirrups w/Fins
#800-0342-R Great White® Robotic Stirrups w/o Fins

Required Accessory: Purchase Two SpringLoc Clamp

#800-0338 **#800-0338-EU**
#800-0338-JPN **#800-0338-SWISS**

Optional Accessories


#508-1354 Replacement Boot Pads w/o Fins
#508-1415 Replacement Boot Pads w/Fins
#800-0284 Stirrup Wall Rack
#800-0074-R Robotic Stirrups Dolly
#800-0160 Disp. Clear View Leg Drapes,
Sterile, 10/cs*

*Sterile products are non-returnable.



Patent Pending

600 lb. (272 kg)
Patient Weight Capacity



#800-0342-RF

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351