

Surgical Peg Board System

Easily secure your patient to precise positions


Our Peg Board for surgical procedures is a durable single board that does not have to be assembled and broken down between cases.

Includes:

- (4) 6" covered short pegs
 - (4) 9" covered long pegs
 - (2) Fixation clamps
 - Gel Overlay item #508-0440
 - Requires (2) SpringLoc Clamps to secure Peg Board to OR table
 - 48"L x 19"W x 1" (119 cm x 48 cm x 3 cm)
- #800-0354** Surgical Peg Board System

Required Accessory: Purchase Two SpringLoc Clamp

#800-0338
 #800-0338-JPN
 #800-0338-EU
 #800-0338-SWISS



Replacement Parts

- #800-0354-A** 6" (15 cm) Removable Pegs, Plastic, Set of Four
- #800-0354-B** 9" (23 cm) Removable Pegs, Plastic, Set of Four
- #800-0354-C** Fixation Clamps, Set of Two
- #800-0354-D** Peg Board/Posting Board
- #800-0354-E** 6" (15 cm) Removable Pegs, Stainless Steel, Set of Four
- #800-0354-F** 9" (15 cm) Removable Pegs, Stainless Steel, Set of Four
- #508-0440** Gel Overlay, 20" x 46" x .5" (51 cm x 117 cm x 1.3 cm)
- #800-0354-G** Gel Overlay 20" x 46" x .5" (61 cm x 117 cm x 1.3 cm)
- #800-0354-H** 4" (10 cm) Peg Extension, ea

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351