

**NEW!**

# Great White® HD 1000 Stirrup Boot System

Discover innovation in our interchangeable stirrup boot system with 1,000 lb Patient Wt. Capacity

Includes 3  
Different  
Style  
Boots!

## Features

- **Proprietary 3-in-1 System** includes Platinum, Bariatric Maxima & Premium Stirrup Boots
- **Stronger!** A new, engineered stainless steel tube increases load capacity to 1,000 lb
- **Redesigned!** 3-axis control handle mimics the human hip's motion for smooth and accurate high/low lithotomy, adduction & abduction positioning
- **Locking Gas Spring** locks stirrup in place—stays locked until trigger is pulled
- **Ranges** Lithotomy Range -25° low to 75° high; Adduction -9° and abduction +25°
- **Included** Premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- **Reduced Setup** Quick, easy connection reduces setup time, training and downtime
- **Easily Repositions** Lithotomy indicators accurately reposition legs

#800-0364 Great White® HD 1000 Stirrup Boot System

**Required Accessory: Purchase Two SpringLoc Clamps**

#800-0338

#800-0338-EU

#800-0338-JPN

#800-0338-SWISS

## Optional Accessory

#800-0160 Disp. Clear View Leg Drapes, Sterile, 10/cs\*

\*Sterile products are non-returnable.

## Replacement Pads

#508-1415 Replacement Boot Pads w/Fins

#508-1354 Replacement Boot Pads w/o Fins

#508-1502 Replacement Bariatric Maxima Boot Pads



Select from the included (A) Bariatric Maxima, (B) Platinum & (C) Premium Boots.



Quick and easy connection reduces setup time making it a cost effective product.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351