



# Ortho BioMetrix (KRS) Knee Arthroscopy System

## *Innovative advancement in patient positioning*

The unique lightweight design conforms to patients legs of all sizes, while providing safe and secure positioning that allows for full access of the knee joint from all angles. Inflatable and disposable air bladders assure conformity to any patient leg.

### **Features**

- *Easy to set-up and use*
- *Total access to knee joint, allowing stress & manipulation of joint during procedure*
- *Increased Access allows surgeon to perform surgery while sitting, lessening fatigue*
- *Extremely lightweight, weighing only 1.5 lbs*
- *Enhanced efficacy, safety & comfort*
- *Facilitates sterile draping*
- *Promotes adjust-ability & control*
- *Provides support for Distal Thigh*
- *Simple to remove for additional procedures*
- *Made from polycarbonate materials, the KRS is a strong and long lasting device that can be cleaned with any hospital grade disinfectant*



- *If High Level Disinfection is desired, the KRS can be disinfected with an approved Glutaraldehyde or similar solution*
- *Sterilization can also be achieved through the standard Autoclave Process using steam sterilization*

### **#810-0001**

Ortho BioMetrix Knee Arthroscopy System

### **Optional Accessories**

**#810-0002** Repl. Clamps & Brackets, 10/cs

**#810-0003** Repl. Air Bladders, 25/cs

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351