


Coated Positioning Pads

Features

- Reusable antimicrobial coated positioners ideal for extremity positioning
- The positioners are radiolucent and the coating is impregnated with Agion antimicrobial ionic silver that inhibits the growth of bacteria
- Reduces chance of cross-contamination
- Slip-Resistant



Item#	Description	Dimensions
A. #508-0355	Double Angle Wedge, 35°/55°, ea	24"L x 9.75"W x 6.75"H (61cm x 25cm x 17cm)
#508-0356	Angle Wedge, 45°, ea	24"L x 7"W x 7"H (61cm x 18cm x 18cm)
#508-0357	Wedge, 15°, ea	10.5"L x 9.75"W x 2.75"H (27cm x 25cm x 7cm)
B. #508-0358	Double Angle, 55°/35°, ea	14"L x 9.5"W x 6.5"H (36cm x 24cm x 17cm)
#508-0359	Wedge, 20°, ea	11"L x 8.5"W x 3.75"H (28cm x 22cm x 10cm)
#508-0360	Head & Neck, ea	10.75"L x 7"W x 6"H (27cm x 18cm x 15cm)
C. #508-0361	Head Support, Occipital Adult, ea	10"L x 10"W x 4"H (25cm x 25cm x 10cm)
D. #508-0364	Semicircle, ea	14"L x 8.5"W x 4.5"H (36cm x 22cm x 11cm)
#508-0365	Roll, Small, ea	24"L x 6"D (61cm x 15cm)
#508-0366	Roll, Large, ea	24"L x 8"D (61cm x 20cm)
E. #508-0367	Donuts, Large, 6 pk.	9"D x 1.5"H (23cm x 4cm)
#508-0368	Donuts, Small, 6 pk.	6"D x 1.5"H (15cm x 4cm)
F. #508-0362	Rectangular, Large, ea	14"L x 9.25"W x 2"H (36cm x 24cm x 5cm)
#508-0363	Rectangular, Small, ea	9.5"L x 7"W x 2"H (24cm x 18cm x 5cm)

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | www.schuremed.com

orders@schuremed.com | 452 Randolph Street, Abington, MA 02351