



## Schure Slide XLA Ambulance

*Choose Flexible, Rigid or Full Size Rigid*

The Schure Slide XLA Ambulance is an ultra-smooth flexible sheet used to transfer patients safely and easily from bed to stretcher and/or from stretcher to bed. The sliding action eliminates the need to lift the patient, removing the risk of possible back injury to emergency personnel.

### **Features**

- *Unique design allows patient transfer in a sitting or supine position in comfort*
- *Ergonomically-designed handles enhance slide placement & removal*
- *Non-disposable, radiolucent, stain-resistant – cleans with hospital-approved disinfectant*
- *Stores conveniently under stretcher mattress for quick accessibility*
- *Works effectively with cotton draw sheets or disposable sheets*
- *Potential savings on Workman's Compensation Insurance*

**#800-0097** Flexible, 33"L x 18"W x 1/16"

**#800-0166** Rigid, 33"L x 18"W x 3/16"

**#800-0126** Full-Size Rigid, 72"L x 18"W x 3/16"



**SCHURE SLIDE**  
PATIENT TRANSFER SYSTEM

**Important!** Please note that the Schure Slides are made for sliding a patient from one surface to another and are not constructed to support patient weight.

Phone 888.724.8763 • 781.982.7000 | Fax 781.982.7001 | [www.schuremed.com](http://www.schuremed.com)

[orders@schuremed.com](mailto:orders@schuremed.com) | 452 Randolph Street, Abington, MA 02351