

E-Z LIFT BEACH CHAIR



#800-0142

500 lb (227 kg)
Patient Weight Capacity

Our double-ball joint head positioning system safely secures the patient's head even with kyphosis patients



Shoulder pull-away panels provide excellent surgical exposure

A light squeeze of the trigger engages the powerful piston system to safely lift or lower the patient



Any surgical OR table can be converted to a Beach Chair for arthroscopic surgery



E-Z LIFT BEACH CHAIR



Position patients for shoulder surgery with the lightweight shoulder arthroscopy beach chair from SchureMed

Position patients for shoulder surgery with this lightweight and cost-effective surgical beach chair that adjusts from 0-90° with a powerful piston system that assists O.R. staff when lifting the patient into position. A safe and secure option to our fully-powered shoulder surgery chair system.

Features

- ◆ Significantly reduces effort by O.R. staff when lifting patients over other manual systems
- ◆ Double-ball joint head system safely secures patient's head
- ◆ Sliding head rest conforms to every patient
- ◆ Articulation range 0-90°
- ◆ 500 lb. surgical patient weight capacity
- ◆ 28"W x 30"D x 47"H, 32 lbs
- ◆ Five-piece Deluxe Pad Set & Lateral Brace included

ACCESSORIES

For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at orders@schuremed.com or visit us online at www.schuremed.com



#800-0142 E-Z Lift Beach Chair

Required Accessory:

Select One Model

#800-0164 Disp. Face Masks, 12/cs

#800-0369 Disp. Head & Chin Straps, 12/cs

**All accessories listed for the E-Z Lift Beach Chair will also work with the Powered Beach chair*

Optional Accessories

#800-0074 Shoulder Chair Dolly

#800-0050 Multi-Axis Arm Positioner

#800-0134 Schure Socket XL

#800-0294 Complete Shoulder System w/
Schure Loc XPS & E-Z Lift Beach Chair