

HD 1000 STIRRUPS

#800-0364



Great White
HD 1000
Stirrup Boot System

1,000 lb (454 kg)
Patient Weight Capacity

Lightly squeeze the ergonomic handle to easily allow full range of interoperative patient positioning

Simply pull plungers to securely lock boots into place

Locking gas spring locks stirrup in place - and stays locked until the trigger is pulled. This allows for superior patient positioning for obese patients

3-Axis control

Required Accessory
#800-0338
Available in EU-JPN-SWISS



For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at orders@schuremed.com or visit us online at www.schuremed.com

HD 1000 STIRRUPS



A Bariatric Maxima #508-1502

B Platinum #508-1415

C Premium #508-1354



Quick and easy connection reduces setup time, training, and downtime.

Features

- ◆ Our proprietary 3-in-1 System includes three pairs of Stirrup Boots: Premium, Platinum & Bariatric
- ◆ Locking Gas Spring locks stirrup in place - stays locked until trigger is pulled
- ◆ Lithotomy Range -25° low to 75° high; Adduction 9° and abduction 25°
- ◆ Quick, easy connection reduces setup time, training and downtime
- ◆ Premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning included.
- ◆ 3-Axis control handle mimics the human hip's motion for smooth and accurate high/low lithotomy, adduction, and abduction positioning
- ◆ The included premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- ◆ Engineered stainless steel tube increases patient weight capacity to 1,000lbs (454kg)

#800-0364 Great White HD 1000 Stirrup Boot System

Required Accessory:

Purchase Two SpringLoc Clamps

#800-0338 #800-0338-JPN
#800-0338-EU #800-0338-SWISS

Optional Accessories

#800-0160 Disp. Clear View Leg Drapes

Replacement Pads

#508-1415 Replacement Boot Pads w/ Fins
#508-1354 Replacement Boot Pads w/o Fins
#508-1502 Repl. Bariatric Maxima Boot Pads