

# SCHURESOCKET INFINITY

\*patent pending



#800-0386  
#800-0386-SQ


#800-0386-F

www.schuremed.com

## THE INFINITY

Patient weight capacity is 800 lbs (363 kg)

Accepts:  
5/8 Rod  
.709 (18mm) Rod  
3/4 Rod  
1/4" x 1" Flat Bar  
3/8" x 1" Flat Bar




The Infinity Line's Adjustable Jaw Allows The Infinity Clamps To Attach To Any OR Table In The World!

Eliminate the need to purchase different clamps for different surgical table side rail sizes with this innovative new design.

## INFINITYsq

Patient weight capacity is 800 lbs (363 kg)


Designed to accept .787" (20mm) Square Bar



## INFINITY BLADE CLAMP

Patient weight capacity is 800 lbs (363 kg)

Accepts:  
Great White Stirrup Blade  
5/8 Rod  
.709 (18mm) Rod  
3/4 Rod  
1/4" x 1" Flat Bar  
3/8" x 1" Flat Bar



Quickly and easily connects and disconnects securely anywhere along any side rail, even over side rail stand off.

For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at [orders@schuremed.com](mailto:orders@schuremed.com)