

PEDIATRIC STIRRUPS



#800-0342-PUPS | #800-0342-KIDS

160 lb (73 kg)
PUPS Weight Capacity

250 lb (113 kg)
KIDS Weight Capacity

Lightly squeeze the ergonomic handle to easily allow full range of interoperative patient positioning



Required
Accessory
#800-0338
Available in
EU-JPN-
SWISS



PEDIATRIC STIRRUPS

Comfortable and secure leg positioning for pediatric patients during urological and laparoscopic procedures.

The Pediatric E-Z Lift Stirrups are as easy-to-use as its forerunner. Exclusive axial rotation from the handle for superior patient positioning. Boot adjusts with a single clamp which simplifies patient positioning. Boot pads included.

Features

- ◆ E-Z Lift feature easily control the patient's leg when moving it into the operative position
- ◆ Floating boot automatically adjusts when repositioning the pediatric boot stirrup
- ◆ Range: Abduction +25°, Adduction -9°, Lithotomy range of +75° to -30° & externally rotates out 20°
- ◆ Exclusive axial rotation from the handle for superior patient positioning
- ◆ #800-0342-PUPS 160 lb (73 kg) patient capacity
- ◆ #800-0342-KIDS 250 lb (113 kg) patient capacity

ACCESSORIES



#800-0342-PUPS & #800-0342-KIDS Pediatric Stirrups
#800-0238 Pediatric E-Z Lift Stirrups Bundle (Includes 1 pair of each stirrup)

Required Accessory:

Purchase Two SpringLoc Clamps

#800-0338 #800-0338-JPN
#800-0338-EU #800-0338-SWISS

Replacement Pads

#508-1295 Replacement PUPS Boot Pads
#508-1296 Replacement KIDS Boot Pads

Optional Accessories

#800-0074-P Pediatric Stirrup Dolly
#800-0160 Disp. Clear View Leg Drapes

For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at orders@schuremed.com or visit us online at www.schuremed.com