

PLATINUM STIRRUPS



#800-0342-PL

600 lb (272 kg)
Patient Weight Capacity

Lightly squeeze the ergonomic handle to easily allow *full range of interoperative patient positioning*



Locking gas spring locks stirrup in place - and stays locked until the trigger is pulled. This allows for superior patient positioning for obese patients

3-Axis control



Required Accessory #800-0338 Available in EU-JPN-SWISS

PLATINUM STIRRUPS



Gynecology, laparoscopy & urology procedures are safer than ever

The Great White Platinum Stirrups provide safe and quick positioning assistance for improved surgical site access. This intra-operative, adjustable medical stirrup for gynecological, laparoscopic, and urological applications offers exceptionally precise placement, effectively minimizing the possibility of nerve damage. Exclusive axial rotation, high/low lithotomy and abduction/adduction from the handle improves patient positioning.

Features

- ◆ Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- ◆ Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- ◆ Lithotomy Range -30° low to 75° high; Adduction 9° and abduction 25°
- ◆ The raised lateral fin boots effectively minimize the possibility of nerve damage
- ◆ Included Premier boot liners wrap around the foot, ankle and calf providing safe, comfortable and secure patient positioning
- ◆ Lithotomy indicators help to accurately reposition legs

ACCESSORIES



Wall
Rack
#800-0284



Stirrup
Dolly
#800-0074-S



Replacement
Pads
#508-1415

#800-0342-PL Great White Platinum Stirrups

Required Accessory:

Purchase Two SpringLoc Clamps

#800-0338 #800-0338-JPN
#800-0338-EU #800-0338-SWISS

Replacement Pads

#508-1415 Replacement Boot Pads w/ Fins
#508-1354 Replacement Boot Pads w/o Fins

Optional Accessories

#800-0074-S Great White Stirrup Dolly
#800-0284 Great White Stirrup Wall Rack
#800-0160 Disp. Clear View Leg Drapes

For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at orders@schuremed.com or visit us online at www.schuremed.com