

# POWERED BEACH CHAIR

#800-0004



**500 lb (227 kg)**  
Patient Weight Capacity

Our double-ball joint head positioning system safely secures the patient's head even with kyphosis patients

A self-contained motor independent of surgical table power automatically lifts patients - accommodates patients up to 500lbs (227kg)

Any surgical OR table can be converted to a Beach Chair for arthroscopic surgery

Shoulder pull-away panels provide excellent surgical exposure

90 to 0°



# POWERED BEACH CHAIR



#800-0004 Powered Beach Chair

## ***The only powered surgical beach chair on the market with a self-contained motor***

The SchureMed Powered Beach Chair is the only shoulder surgery positioning chair on the market with a self-contained motor that does not rely on the surgical table for power. Our double-ball joint head positioning system safely secures the patient's head even with kyphosis patients. Attaches to ANY operating room table with our integrated, adjustable mounting clamps. Lateral brace included.

### **Features**

- ◆ Keeps the patient and operating room staff safer and more secure
- ◆ Attaches to ANY operating room table with our integrated, adjustable mounting clamps
- ◆ 5-Piece Deluxe Pad Set and Lateral Brace included
- ◆ 1,350 lbs. of lift securely supports arthroscopy patients up to 500 lbs
- ◆ Shoulder pull-away panels provide excellent surgical exposure
- ◆ Articulation range 0-90° in 8 seconds in a slow, controlled movement

## **COMPLETE SHOULDER SYSTEM**

Complete Shoulder System includes the Schure Loc™ (XPS #800-0280), Premium Head & Chin Straps (#800-0369), Multi-Axis Arm Positioner (#800-0050), Schure Socket XL (#800-0134), Disposable Schure Loc™ XPS Sterile Kit (#800-0281), the Shoulder Chair Dolly (#800-0074), and either the Powered Beach Chair (#800-0004) or the E-Z Lift Beach Chair (#800-0142).



For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at [orders@schuremed.com](mailto:orders@schuremed.com) or visit us online at [www.schuremed.com](http://www.schuremed.com)