

PREMIUM STIRRUPS

#800-0342-PR



400 lb (181 kg)
Patient Weight Capacity

Lightly squeeze the ergonomic handle to easily allow full range of interoperative patient positioning



Required
Accessory
#800-0338
Available in
EU-JPN-
SWISS



PREMIUM STIRRUPS



High quality, feature-rich, and priced right, the Great White Premium Stirrups are exactly what your surgical team requires. The enhanced range of motion allows for -30° for low lithotomy, to 75° high lithotomy; adduction 9° to abduction 25° and also externally rotates out 20°. Two required SpringLoc Clamps sold separately.


Features

- ◆ Exclusive 3-axis control from the handle mimics the anatomy of the human hip
- ◆ Simply squeeze the trigger for smooth and accurate high/low lithotomy, adduction and abduction positioning
- ◆ Lithotomy Range -30° low to 75° high; Adduction 9° and abduction 25°
- ◆ Included premier boot liners wrap the foot, ankle & calf providing safe, comfortable & secure patient positioning
- ◆ Lithotomy indicators help to accurately reposition legs
- ◆ Lift assist allows easy control with minimal effort
- ◆ Flexible boot accommodates larger patients

ACCESSORIES



Wall
Rack
#800-0284



Stirrup
Dolly
#800-0074-S



Replacement
Pads
#508-1415

#800-0342-PR Great White Premium Stirrups

Required Accessory:

Purchase Two SpringLoc Clamps

#800-0338 #800-0338-JPN
#800-0338-EU #800-0338-SWISS

Replacement Pads

#508-1415 Replacement Boot Pads w/ Fins
#508-1354 Replacement Boot Pads w/o Fins

Optional Accessories

#800-0074-S Great White Stirrup Dolly
#800-0284 Great White Stirrup Wall Rack
#800-0160 Disp. Clear View Leg Drapes

For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at orders@schuremed.com or visit us online at www.schuremed.com