



SchureMed™ Surgical Stirrup Comparison Chart

	PRODUCT	Patient Weight Limit	Protective Fin Design	Lift Mechanism	Lithotomy Adjustment (From Sterile Field)	Lithotomy Range (From Sterile Field)	Abduction Adjustment (From Sterile Field)	Abduction Range (From Sterile Field)	Lithotomy Indicator
BEST	Great White® HD 1000 Stirrups #800-0364	1,000 lbs. (454 kg)	✓	✓	Yes, All Positions In Range	+75°, -25° External Rotation 20°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
	Great White® Maxima Stirrups #800-0342-M	800 lbs. (363 kg)	✓	✓	Yes, All Positions In Range	+75°, -30° External Rotation 20°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
	Great White® Platinum Stirrups #800-0342-PL	600 lbs. (272 kg)	✓	✓	Yes, All Positions In Range	+75°, -30° External Rotation 20°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
	Great White® Robotic Stirrups #800-0342-RF	600 lbs. (272 kg)	✓	✓	Yes, All Positions In Range	+60°, -55° External Rotation 20°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
BETTER	Great White® Premium Stirrups #800-0342-PR	400 lbs. (181 kg)	X	✓	Yes, All Positions In Range	+75°, -30° External Rotation 20°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
GOOD	Lithotomy Stirrups #800-0049	350 lbs. (159 kg)	✓	X	X	N/A	N/A	N/A	N/A
	Candy Cane Stirrups #800-0012	350 lbs. (159 kg)	✓	X	X	N/A	N/A	N/A	N/A
PEDIATRIC	Pediatric E-Z Lift Stirrups #800-0342-PUPS #800-0342-KIDS	PUPS: 160 lbs. (73 kg) KIDS: 250 lbs. (113 kg)	X	✓	Yes, All Positions In Range	+75°, -30° External Rotation 20°	Yes, All Positions In Range	+25°, -9°	✓
	Pediatric Lithotomy Stirrups #800-0233-PUPS #800-0234-KIDS	PUPS: 160 lbs. (73 kg) KIDS: 250 lbs. (113 kg)	X	X	X	N/A	N/A	N/A	N/A