



# TKR POSITIONER FOR TOTAL KNEE REPLACEMENT



The ball track system at the heel of the leg holder places the operative leg in flexion, extension, tilt, and rotation.



Smooth sliding motion for effortless leg positioning



The required CamLoc Clamp is available for UK, EU, DEN, SWISS, and JPN rails



Steam Sterilization Tray Included!

# TKR POSITIONER FOR TOTAL KNEE REPLACEMENT



**Highly functional, easy to adapt positioner for total knee replacement surgery. Completely secure and dynamic range of leg positions.**

This ball-track system for total knee replacements is a huge leap forward for the surgeon and patient alike. Ball joint at the heel of the leg holder places the operative leg in extension, flexion, tilt and rotation. Easy-to-use and completely rigid once locked. Eliminates the need for additional leg positioning devices, sandbags and gel-pads.

## Features

- ◆ Custom Sterilization Tray and 1 TKR Boot Liner included with purchase
- ◆ Track Locking Bracket securely locks the Base Plate to the table with two thumb screws
- ◆ Locks in any position with a simple turn of the wrist
- ◆ Smooth sliding motion for effortless leg positioning
- ◆ Well leg lies flat – eliminating the need for well-leg positioners
- ◆ All components are autoclavable (Exposure Time at 270°F (132°C) for 4 minutes)
- ◆ Available in 24" (61 cm) & 30" (76 cm)

## ACCESSORIES



**#800-0141 TKR Positioner for Total Knee Replacement, 24" (61cm)**  
**#800-0171 TKR Positioner for Total Knee Replacement, 30" (76cm)**

### Required Accessory:

#### Select One CamLoc Clamp

#800-0271 for use with 1-2 drapes  
#800-0271-L for use with 3-5 drapes  
#800-0271-UK #800-0271-EU  
#800-0271-DEN #800-0271-JPN  
#800-0271-SWISS

### Optional Accessories

#800-0146 Disposable TKR Boot Liners,  
Sterile, 12/cs

For more information, or to place an order, contact our friendly and knowledgeable customer service team at +781-982-7000 or email us at [orders@schuremed.com](mailto:orders@schuremed.com) or visit us online at [www.schuremed.com](http://www.schuremed.com)