

## Stirrup Clamp

*Instantly attaches to side rail quickly & smoothly*

The newly designed Stirrup Clamp locks patient positioning devices and surgical equipment with a flat mounting post anywhere along an operating table side rail.

### **Features**

- Securely locks accessories with a flat 1" x ¼" (2.5cm x .6 cm) mounting post
- Attaches easily to the operating table side rail without the need of table notches
- Stirrup Clamp is easily removed without the need for a release lever
- 1,000 lb. (454 kg) patient weight capacity



### **Stirrup Clamp**

**#800-0396**

**#800-0396-EU**

*Each clamp sold separately.*