

Surgical Beach Chairs

#800-0004 & #800-0142

Our double-ball joint head positioning system safely secures the patient's head even with kyphosis patients

A light squeeze of the trigger engages the powerful piston system to safely lift or lower the patient

Any surgical OR table can be converted to a Beach Chair for arthroscopic surgery

**500 lb. (227 kg)
Patient Weight
Capacity**

#800-0142
E-Z Lift
Beach Chair

Shoulder pull-away panels provide optimal surgical exposure

#800-0004
Powered
Beach Chair

A self-contained motor independent of surgical table power automatically lifts patients - accommodates patients up to 500lbs (227kg)



International Part Numbers:

#800-0142

#800-0142-DEN

#800-0142-JPN

#800-0142-EU

#800-0142-UK

#800-0142-SWISS

INCLUDED ACCESSORIES



#508-1461
Shoulder Chair
Pad Set



#800-0091
Shoulder Chair
Lateral Brace



#790-0063
Shoulder Chair
Head Section



#800-0369
Premium
Head & Chin Straps

Included:

1x Set of Deluxe Shoulder Chair Pads

1x Shoulder Chair Lateral Brace

1x Shoulder Chair Head Section

1x Premium Head & Chin Straps

OPTIONAL ACCESSORIES



#800-0074
Shoulder Chair
Dolly



#800-0398
The Ótimo
Armboard



#800-0164
Disposable
Face Masks



#790-0103
SchureLoc XPS
Retractor Holder

#790-0040
SchureLoc XPS
Forearm
Attachment

#790-0039
SchureLoc XPS
Adapter



#800-0281
Disposable
SchureLoc XPS
Sterile Kit



#800-0280
SchureLoc XPS

Optional Accessories

#800-0074 Shoulder Chair Dolly

#800-0398 The Ótimo Armboard

#800-0164 Disposable Face Masks

#800-0280 SchureLoc XPS

#790-0039 SchureLoc XPS Adapter

#790-0040 SchureLoc XPS Forearm Attachment

#800-0281 SchureLoc XPS Sterile Kit