



Great White
HD 1000
Stirrup Boot System



**Three Stirrups in One.
Every Patient Covered.**



Stainless steel tube
for increased patient
weight capacity of
1,000lbs

Pull the plunger
and lift the boot
to remove



**3-in-1 boot system
with quick connect
improves turnaround times**

Simply drop the boot
in and listen for the
snap for audible
confirmation of
connection



INCLUDED ACCESSORIES

Included:

1x Set Premium Boots w/ pads

1x Set of Platinum Boots w/ Pads

1x Set of Maxima Boots w/ pads

#800-0364-PR
Premium
Boots w/ Pads



#800-0364-PL
Platinum
Boots w/ Pads



#800-0364-M
Maxima Boots
w/ Pads



OPTIONAL ACCESSORIES

Required Accessories: **2x Stirrup Clamps**

#800-0396
Stirrup Clamp



#800-0074-HD1000
HD1000 Dolly



#800-0284
Stirrup Wall Rack



#508-1354
Premium Boot
Pads



#508-1502
Bariatric Boot
Pads



#508-1415
Platinum Boot
Pads



Optional Accessories:

#800-0074-HD1000 HD1000 Stirrup Dolly

#800-0284 Stirrup Wall Rack

#508-1354 Premium Boot Pads

#508-1415 Platinum Boot Pads

#508-1502 Maxima Boot Pads