

Maxima Stirrups

#800-0342-M



Boot free-floats throughout each range reducing the risk of superficial nerve injury

Axial rotation from the handle for superior patient positioning

A wide boot and a stronger gas cylinder make this stirrup perfect for your heaviest patients

Extended lateral boot design helps reduce the risk of peroneal nerve injury

Included:
1x Set Bariatric Boot Pads

#800-0074-S
Stirrup Dolly

#508-1502
Bariatric Boot
Pads

Optional Accessories:
#800-0074-S Stirrup Dolly
#800-0284 Stirrup Wall Rack
#508-1502 Bariatric Boot Pads

Required Accessories:
2x Stirrup Clamps

#800-0396
Stirrup Clamp



#800-0284
Stirrup Wall Rack

