

Platinum Stirrups

#800-0342-PL

Exclusive axial control
from the handle
mimics the anatomy
of the human hip

Simply squeeze the trigger for
smooth and accurate
high/low lithotomy, adduction
and abduction positioning

The raised lateral fin
boots effectively
minimize the possibility
of nerve damage

Lithotomy indicators
help to accurately
reposition legs

Included:
1x Set Platinum Boot Pads

Required Accessories:
2x Stirrup Clamps

#800-0396
Stirrup Clamp



#508-1415
Platinum Boot
Pads



Optional Accessories

#800-0074-S Stirrup Dolly

#800-0284 Stirrup Wall Rack

#508-1415 Platinum Boot Pads

#800-0284
Stirrup Wall Rack



#800-0074-S
Stirrup Dolly

