

SCHURELOCKXPS

#800-0280

**Newest evolution is
lower profile, lighter,
and provides more
rigid fixation**



#800-0398

The Ótimo enables secure
holding of the patient's arm
during procedure



Adapter joint fits
different limb holders
for superior patient
positioning



Single lever releases
and relocks the
device for safe and
easy repositioning



Clamp fits ALL global
OR table side rails



Six joints for unlimited
intra-operative free re-positioning
in all 3 planes



INCLUDED ACCESSORIES



Included:

1x SchureLoc XPS Adapter
1x XPS Forearm Attachment
1x Sterile Drapes
1x Foam Hand/Forearm Wrap
1x Roll of Self-Adherent Coban

OPTIONAL ACCESSORIES



Optional Accessories

#790-0056 SchureLoc XPS Armboard
 #790-0057 SchureLoc XPS Counter Traction Post
 #790-0075 SchureLoc XPS Elbow/ Ankle Positioner
 #790-0077 SchureLoc XPS TKR Positioner

New Attachments

#790-0101 SchureLoc XPS Femoral Hook
 #790-0103 SchureLoc XPS Retractor Holder
 #790-0058 SchureLoc XPS Finger Trap

Disposables

#800-0281 10x SchureLoc XPS Sterile Kit
 #800-0146 12x Disposable TKR Boot Liners