



We put you and your patients in a better position



STIRRUP SOLUTIONS

888-724-8763 781-982-7000 | www.schuremed.com
452 Randolph Street, Abington, MA 02351



Established in 1994, SchureMed was founded by company President and CEO, Peter Schuerch. Peter's talent for problem solving led him from the operating room to the board room. Peter is not only the President and CEO, but also an integral part of the design team.

SchureMed now manufactures high quality surgical assistance devices for medical facilities all over the world. Up or down, big or small, SchureMed is putting you, and your patients, in a better position.



OUR MISSION

Our mission is to enhance the operating room experience of patients and medical professionals worldwide. Our passion for innovation, quality, and safety drive us to create highly functional surgical positioning solutions that optimize results and reduce operating time through greater efficiency.

OUR VISION

Our vision is to improve quality of life for all patients and medical professionals by becoming the world's most trusted leader in innovative surgical positioning.



TABLE OF CONTENTS



Comparison Chart	02
04	SchureClean
HD 1000 Stirrup System	05
07	Platinum Stirrups
Robotic Stirrups	08
10	Split Leg Stirrups
Maxima Stirrups	11
12	Premium Stirrups
Pediatric Stirrups	13
14	Traditional Stirrup Solutions

SchureMed™ Surgical Stirrup Comparison Chart

	PRODUCT	Patient Weight Limit	Protective Fin Design	Lift Mechanism	Lithotomy Adjustment (From Sterile Field)	Lithotomy Range (From Sterile Field)	Abduction Adjustment (From Sterile Field)	Abduction Range (From Sterile Field)	Lithotomy Indicator
A	Great White® HD 1000 Stirrups #800-0364	1,000 lbs (454 kg)	✓	✓	Yes, All Positions In Range	+75°, -25°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
B	Great White® Maxima Stirrups #800-0342-M	800 lbs (363 kg)	✓	✓	Yes, All Positions In Range	+75°, -30°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
C	Great White® Platinum Stirrups #800-0342-PL	600 lbs (272 kg)	✓	✓	Yes, All Positions In Range	+75°, -30°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
D	Great White® Robotic Stirrups #800-0342-RF	600 lbs (272 kg)	✓	✓	Yes, All Positions In Range	+60°, -55°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
E	Great White® Premium Stirrups #800-0342-PR	400 lbs (181 kg)	X	✓	Yes, All Positions In Range	+75°, -30°	Yes, All Positions In Range	+25°, -9°	✓ (on rod)
F	Lithotomy Stirrups #800-0049	350 lbs (159 kg)	✓	X	X	N/A	N/A	N/A	N/A
G	Candy Cane Stirrups #800-0012	350 lbs (159 kg)	✓	X	X	N/A	N/A	N/A	N/A
H	Pediatric E-Z Lift Stirrups #800-0342-PUPS #800-0342-KIDS	PUPS: 160 lbs (73 kg) KIDS: 250 lbs (113 kg)	X	✓	Yes, All Positions In Range	+75°, -30°	Yes, All Positions In Range	+25°, -9°	✓
I	Pediatric Lithotomy Stirrups #800-0233-PUPS #800-0234-KIDS	PUPS: 160 lbs (73 kg) KIDS: 250 lbs (113 kg)	X	X	X	N/A	N/A	N/A	N/A

GREAT WHITE STIRRUPS



STIRRUPS



PEDIATRIC






Great White HD 1000 Stirrup Boot System



Three Stirrups in One Every Patient Covered



Stainless steel tube
for increased patient
weight capacity of
1,000 lbs (454 kg)

Pull the plunger
and lift the boot
to remove



**3-in-1 boot system
with quick connect
improves turnaround times**

Simply drop the boot
in and listen for the
snap for audible
confirmation of
connection



INCLUDED ACCESSORIES

Included:

1x Set Premium Boots w/ Pads

1x Set of Platinum Boots w/ Pads

1x Set of Maxima Boots w/ Pads

#800-0364-PR
Premium
Boots w/ Pads



#800-0364-PL
Platinum
Boots w/ Pads



#800-0364-M
Maxima Boots
w/ Pads



OPTIONAL ACCESSORIES

Required Accessories:

800-0396 2x Stirrup Clamps

#800-0396
Stirrup Clamp



#508-1354
Premium Boot
Pads



#508-1415
Platinum Boot
Pads



#800-0074-HD1000
HD1000 Dolly



#800-0284
Stirrup Wall Rack



#800-0364-PL-SC
SchureClean
Platinum Boots w/
Pads



#508-1502
Bariatric Boot
Pads



Optional Accessories:

#800-0074-HD1000 HD 1000 Stirrup Dolly

#800-0284 Stirrup Wall Rack

#508-1354 Premium Boot Pads

#508-1415 Platinum Boot Pads

#508-1502 Bariatric Boot Pads

#800-0364-PL-SC SchureClean Platinum Boots w/ Pads

NEW OPTIONAL ACCESSORY
FOR THE HD 1000 STIRRUP

SchureClean Platinum Boots

#800-0364-PL-SC

**SchureClean Incorporates a New Fabric-Free Boot
Designed for Superior Infection Control**



Enhanced Safety:

Engineered with no exposed hard plastic edges and a lateral fin to protect patients and improve stability

Patient-Centric Design:

A larger, padded boot provides enhanced comfort and accommodates diverse patient anatomies

Eliminates Hook & Loop Straps:

Incorporates medical-grade silicone straps for secure positioning and ease of cleaning, minimizing infection risk

Robust and Flexible:

Constructed with a strong, non-tear-propagating material for long-term durability and reliability in the OR



Platinum Stirrups

#800-0342-PL

Exclusive axial control
from the handle
mimics the anatomy
of the human hip

Simply squeeze the trigger for
smooth and accurate
high/low lithotomy, adduction
and abduction positioning

The raised lateral fin
boots effectively
minimize the possibility
of nerve damage

Lithotomy indicators
help to accurately
reposition legs

Included:
1x Set Platinum Boot Pads

Required Accessories:
#800-0396 2x Stirrup Clamps

#800-0396
Stirrup Clamp



#508-1415
Platinum Boot
Pads



Optional Accessories:
#800-0074-S Stirrup Dolly
#800-0284 Stirrup Wall Rack
#508-1415 Platinum Boot Pads

#800-0284
Stirrup Wall Rack



#800-0074-S
Stirrup Dolly




Robotic Stirrups

#800-0342-RF (w/ fins)

#800-0342-R (w/o fins)



High lithotomy of up to 60° allows for better access to surgical sites



Designed in collaboration with robotic surgeons to maximize surgical outcomes



INCLUDED ACCESSORIES

Included:

1x Set of Platinum Boots w/ Pads (w/ fins 800-0342-RF)

-OR-

1x Set Premium Boots w/ Pads (w/o fins 800-0342-R)

#508-1415
Platinum Boot
Pads



#508-1354
Premium Boot
Pads



OPTIONAL ACCESSORIES

Required Accessories:

#800-0396 2x Stirrup Clamps

#800-0396
Stirrup Clamp



#800-0074-R
Robotic Stirrup Dolly



#800-0284
Stirrup Wall Rack



#508-1354
Premium Boot
Pads



#508-1415
Platinum Boot
Pads



Optional Accessories:

#800-0074-R Robotic Stirrup Dolly

#800-0284 Stirrup Wall Rack

#508-1354 Premium Boot Pads

#508-1415 Platinum Boot Pads

Split Leg Stirrups

#800-0342-SL



Easily attaches to any
OR table side rail

Fixed side rail offers
more space for
additional accessories

Lightweight radiolucent
design reduces OR
fatigue whilst ensuring
precise imaging




High and low lithotomy
can be intraoperatively
adjusted and locked into
place safely with the
squeeze of a trigger



Optional Accessories:
**#800-0324 Split Leg
Foot Boards**

#800-0324
Split Leg
Foot Boards



Required Accessories:
**#800-0228-US-SL
2x Split Leg Clamps**



#800-0228-US-SL
Split Leg Clamp

Maxima Stirrups

#800-0342-M



Boot free-floats throughout each range reducing the risk of superficial nerve injury

Axial rotation from the handle for superior patient positioning

A wide boot and a stronger gas cylinder make this stirrup perfect for your heaviest patients

Extended lateral boot design helps reduce the risk of peroneal nerve injury

Included:
1x Set Bariatric Boot Pads

#800-0074-S
Stirrup Dolly

#508-1502
Bariatric Boot
Pads

Optional Accessories:
#800-0074-S Stirrup Dolly
#800-0284 Stirrup Wall Rack
#508-1502 Bariatric Boot Pads

Required Accessories:
2x Stirrup Clamps

#800-0396
Stirrup Clamp



#800-0284
Stirrup Wall Rack



Premium Stirrups

#800-0342-PRF (w/ fins)
#800-0342-PR (w/o fins)

Simply squeeze the trigger
for smooth and accurate
high/low lithotomy,
adduction and abduction
positioning

Included premier boot
liners wrap the foot, ankle
& calf providing safe,
comfortable & secure
patient positioning

Exclusive axial control from
the handle mimics the
anatomy of the human hip

Lift assist allows easy
control with minimal effort

Required Accessories:

#800-0396 2x Stirrup Clamps



Included:

1x Set of Platinum Boots w/ Pads (w/ fins 800-0342-PRF)
-OR-

1x Set Premium Boots w/ Pads (w/o fins 800-0342-PR)

Optional Accessories:

#800-0074-S Stirrup Dolly

#800-0284 Stirrup Wall Rack

#508-1354 Premium Boot Pads

#508-1415
Platinum Boot Pads
w/ fins

#508-1354
Premium Boot Pads



Pediatric Stirrups

#800-0342-PUPS | #800-0342-KIDS E-Z Lift Stirrups
#800-0233-PUPS | #800-0234-KIDS Lithotomy Stirrups



E-Z Lift Stirrups

#800-0342-PUPS
Patients up to
160 lbs (73 kg)

#800-0342-KIDS
Patients up to
250 lbs (113 kg)

Required Accessories: 2x Stirrup Clamps

#800-0396
Stirrup Clamp



Included:
1x Set PUPS Boot Pads
-OR-
1x Set KIDS Boot Pads

Optional Accessories:

#800-0074-P Stirrup Dolly
#508-1295 PUPS Boot Pads
#508-1296 KIDS Boot Pads



#800-0074-P
Stirrup Dolly



Lithotomy Stirrups

#800-0233-PUPS
Patients up to
160 lbs (73 kg)

#800-0234-KIDS
Patients up to
250 lbs (113 kg)

#800-0006
Schure Socket

Required Accessories:
#800-0006 2x Schure Socket



Traditional Stirrup Solutions

#800-0049 Lithotomy Stirrups
#800-0012 Candy Cane Stirrups

Included:

1x Set Platinum Boot Pads

Optional Accessories:

#800-0074-L Stirrup Dolly

#508-1415 Platinum Boot Pads

#508-1415
Platinum Boot
Pads



#800-0074-L
Stirrup Dolly



Required Accessories:

#800-0006 2x Schure Socket


#800-0006
Schure Socket



Required Accessories:

#800-0134 2x Schure
Socket XL

#800-0134
Schure Socket XL



#508-0453
Gel Candy Cane
Stirrup Sheath



#790-0006
Replacement
Candy Cane
Stirrup Straps



Optional Accessories:

#790-0006 Replacement Straps

#508-0453 Gel Stirrup Sheath

Traditional Stirrup Solutions

#800-0041 Knee Crutches
#800-0096 Uro-Trapper System

Included:

1x Set Knee Crutch Pads



#508-0154
Knee Crutch Pads

#800-0134
Schure Socket XL



Required Accessories:

**#800-0134 2x Schure
Socket XL**



#800-0095
Drain Bags



Optional Accessories:

#800-0095 Drain Bags

#800-0098 Sterile Drain Bags

#800-0098
Sterile Drain Bags



Traditional Stirrup Solutions

#800-0055 Nissen Straps
#800-0118 Bariatric Nissen Straps



Optional Accessory:
#800-0100 1x Shoulder Support

#800-0100
Shoulder Support



Required Accessories:
#800-0228 2x Simple Clamps

#800-0228
Simple Clamps




#800-0055
Nissen Straps



#800-0118
Bariatric Nissen Straps







+1-888-724-8763



orders@schuremed.com



452 Randolph Street
Abington, MA, USA 02351



www.schuremed.com



Rev. 001.24