



CamLoc Rail Clamp

#800-0271-DEN
#800-0271-EU

#800-0271-JPN
#800-0271-UK

#800-0271-SWISS
#800-0271-LP-DEN


#800-0271-LP-EU
#800-0271-LP-JPN

#800-0271-LP-UK
#800-0271-LP-SWISS


Attaches to the operating
table rail at any location

Holds 1-2 drapes between
table side rail and clamp

Accepts accessories with round
posts up to 5/8" (1.6 cm) in
diameter and flat mounting
posts to 1" x 3/8" (2.5 cm x .95 cm)



Available in low profile -
a more sleek and compact design
#800-0271-LP



Available in large - for more drape
capacity (Holds 3-5 drapes between table
side rail and clamp)
#800-0271-L

Designed for use in the sterile field
for when drapes are to be placed
between the clamp and siderail to
maintain sterility; attaching
common accessories to nearly all
OR table designs and brands

For use with:

