


Infinity Line Clamps

#800-0386
#800-0386-F

#800-0386-SQ

Infinity Clamps are Universal

The Infinity Line's adjustable jaw allows Infinity Clamps to attach to any OR table in the world – eliminating the need to purchase different clamps for different surgical table side rail sizes



Infinity

#800-0386

5/8" (16 mm) Diameter Round Rod
.709" (18 mm) Diameter Round Rod
3/4" (19 mm) Diameter Round Rod
3/8" x 1" (.95 cm x 2.5 cm) Square Rod




Infinity sq.

#800-0386-SQ

.787" (20 mm) Square Rod

Infinity Line Clamps have
a 800 lbs (363 kg)
Patient Weight Capacity



#800-0386-F

Infinity Fixed

Stirrup Blade
5/8"D (16 mm) Round Rod
.709" (18 mm) Round Rod
3/4"D (19 mm) Round Rod
1/4" x 1" (.635 cm x 2.5 cm) Flat Bar
3/8" x 1" (.95 cm x 2.5 cm) Flat Bar