



Schure Socket XL

#800-0134
#800-0134-DEN

#800-0134-EU
#800-0134-JPN

#800-0134-UK
#800-0134-SWISS

The revolutionary Schure Socket XL snaps on at any rail location and fits all standard operating room table side rails measuring $\frac{3}{8}$ "W x $1\frac{1}{8}$ "H (.95 cm x 3 cm)

Hardened clamping surface creates maximum clamping force, stopping all movement of operating room accessories

Accepts rectangular mounting posts up to 1" x $\frac{3}{8}$ " (2.5 cm x .95 cm)

Accepts round mounting posts up to $\frac{5}{8}$ " (1.6 cm) in diameter

Patient weight capacity is 800 lbs (363 kg)

Quickly and easily connects and disconnects securely anywhere along the side rail with a quick turn, even over side rail stand off

NEXT
Generation



To connect, place top lip over side rail with logo facing up



Tilt clamp down until bottom lip is underneath side rail and you hear a click



To disconnect, push in on the clamp to unlock



Tilt clamp up and pull away from the side rail

For use with:

