


Schure Spine Frame

#800-0235



Push-button actuation to precisely adjust the width between pads, enabling rapid configuration to accommodate varying patient anatomies and sizes

**500 lbs (227 kg)
Patient Weight
Capacity**

Radiolucent and allows for unrestricted C-arm access

Compatible with any OR table, fits the Jackson Spine Table and Allen Flex Frame



INCLUDED ACCESSORIES



#800-0074-SF
Schure Spine
Frame Dolly



#508-0472
Schure Spine
Frame Pad Set



#800-0277
Schure Spine
Crank Handle



#508-0373
Pair of Disposable
Laminectomy
Arm Cradles



#508-1340
Disposable
Extended Richards
Slotted Headrest

Included:

Crank Handle
Schure Spine Dolly
Schure Spine Frame Pad Set
Disposable Extended Richards
Slotted Headrest
Set of 2 Disposable Arm Cradles


OPTIONAL ACCESSORIES



#800-0074-SF
Schure Spine
Frame Dolly



#508-0472
Schure Spine
Frame Pad Set



#800-0277
Schure Spine
Crank Handle



#508-0373
Case of 12 Pairs of
Disposable
Laminectomy
Arm Cradles



#508-1340
Case of 6
Disposable
Extended
Richards Slotted
Headrests

Replacement Accessories:

#800-0074-SF Schure Spine Dolly
#800-0277 Schure Spine Crank Handle
#508-0472 Schure Spine Frame Pad Set

Disposables:

#508-0373 Case of 12 Pairs of Disposable
Laminectomy Arm Cradles
#508-1340 Case of 6 Disposable Extended
Richards Slotted Headrests