



Simple Clamp

#800-0228
#800-0228-DEN

#800-0228-EU
#800-0228-JPN

#800-0228-UK
#800-0228-SWISS

Designed to secure accessories with a flat mounting blade to most OR table designs and brands

Patient weight capacity is 720 lbs (327 kg)

Securely locks surgical accessories with flat 1" x 1/4" (2.5 cm x .64 cm) horizontally mounting blades

Securely locks surgical accessories with flat 1" x 3/8" (2.5 cm x .95 cm) vertically mounting blades

Slides on to the rail at operating table rail end or notches



For use with:

