

Stirrup Clamp

#800-0396

#800-0396-DEN

#800-0396-EU

#800-0396-JPN

#800-0396-UK

#800-0396-SWISS

Our low-profile Stirrup Clamp is one of the world's easiest-to-use blade clamps, effortlessly securing patient-positioning devices and surgical equipment with a flat mounting post anywhere along an operating-table side rail

Securely
locks accessories
with a flat 1" x 1/4" (2.5
cm x .64 cm)
mounting post

1,000 lbs (454 kg)
patient weight
capacity

Drops effortlessly onto
the operating table side
rail without the need for
table notches

Easily removeable
without the need for a
release lever



For use with:

