




GENOS TKR POSITIONER

#800-0409 30" GENOS TKR Positioner



GENOS TKR POSITIONER

#800-0409 30" GENOS TKR Positioner



The ball track system at the heel of the leg holder allows surgeons to tilt and rotate the operative leg

Pull apart system allows for smaller sterilization tray


Track locking bracket securely locks the Base Plate to the table with two thumb screws

Built-in stops on each end of the track system prevent the leg holder from sliding off of the Base Plate

Locks in any position with a simple turn of the wrist

Place the operative leg in any degree of flexion to full extension

Smooth sliding motion for effortless leg positioning



INCLUDED ACCESSORIES



All components are
autoclavable (Exposure
Time at 270°F (132°C)
for 4 minutes)



#800-0409-TRAY
18" Genos TKR
Sterilization Tray



#800-0146-SGL
Disposable TKR
Boot Liner

Included:
1x TKR Positioner
1x TKR Disposable Boot Liner
1x TKR Sterilization Tray

REQUIRED ACCESSORIES



#800-0271
CamLoc Clamp

Required Accessories
#800-0271 CamLoc Clamp

OPTIONAL ACCESSORIES

Disposables
#800-0146 12x Disposable TKR Boot Liners



#800-0146
Case of 12
Disposable TKR
Boot Liners

Optional Accessories
#800-0271-L CamLoc Clamp Large
#800-0271-LP CamLoc Clamp Low-Profile



#800-0271-L
CamLoc Clamp
Large



#800-0271-LP
CamLoc Clamp
Low-Profile