

Universal Rail Clamp

#800-0085


Designed for a wide-range of purposes; affixing most operating room table accessories to most OR table designs and brands

Accepts any round mounting posts up to 5/8" (1.6 cm) in diameter and rectangular mounting posts up to 1" x 3/8" (2.5 cm x .95 cm)

Patient weight capacity is 350 lbs (159 kg)

Snaps on at any rail location and attaches to all 3/8"W x 1-1/8"H (.95 cm x 3 cm) surgical table side rails

Holds accessories in the vertical and horizontal position



For use with:

